

For Research Use Only

Mouse IL-4 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 85179-10-PBS



Basic Information

Catalog Number:

85179-10-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3750

GenBank Accession Number:

BC027514

GeneID (NCBI):

16189

UNIPROT ID:

P07750

Full Name:

interleukin 4

Purification Method:

Protein A purification

CloneNo.:

251776F7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



EL-4 cells were cultured unstimulated or stimulated with 10 μ g/mL PHA and 10 ng/mL PMA for 2 days. The mean IL-4 concentration was determined to be 49.1 pg/mL in unstimulated EL-4 supernatant, 328.2 pg/mL in stimulated EL-4 supernatant.



Sandwich ELISA standard curve of MP01896-4, Mouse IL-4 Recombinant Matched Antibody Pair - PBS only. 85179-10-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3750. 85179-9-PBS was HRP conjugated as the detection antibody. Range: 15.6-500 pg/mL