

For Research Use Only

SHC4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85181-1-PBS



Basic Information

Catalog Number:

85181-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3336

GenBank Accession Number:

BC033907

GeneID (NCBI):

399694

UNIPROT ID:

Q6S5L8

Full Name:

SHC (Src homology 2 domain containing) family, member 4

Calculated MW:

630 aa, 69 kDa

Purification Method:

Protein A purification

CloneNo.:

242805D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01891-2, SHC4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85181-2-PBS. Detection antibody: 85181-1-PBS. Standard: Ag3336. Range: 0.781-100 ng/mL.