

For Research Use Only

EPS8 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85184-1-PBS



Basic Information

Catalog Number: 85184-1-PBS	GenBank Accession Number: BC030010	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 2059	CloneNo.: 242835B3
Source: Rabbit	UNIPROT ID: Q12929	
Isotype: IgG	Full Name: epidermal growth factor receptor pathway substrate 8	
Immunogen Catalog Number: AG3121	Calculated MW: 822 aa, 92 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01889-2, EPS8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85184-4-PBS. Detection antibody: 85184-1-PBS. Standard: Ag3121. Range: 0.313-40 ng/mL