

For Research Use Only

NDUFA5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85187-3-PBS



Basic Information

Catalog Number:

85187-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9995

GenBank Accession Number:

BC000813

GeneID (NCBI):

4698

UNIPROT ID:

Q16718

Full Name:

NADH dehydrogenase (ubiquinone) 1

alpha subcomplex, 5, 13kDa

Calculated MW:

13 kDa

Purification Method:

Protein A purification

CloneNo.:

242860D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01899-2, NDUFAS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85187-3-PBS. Detection antibody: 85187-2-PBS. Standard: Ag9995. Range: 0.781-100 ng/mL.