

For Research Use Only

# NSMCE2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85188-2-PBS



## Basic Information

Catalog Number:

85188-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4545

GenBank Accession Number:

BC032797

GeneID (NCBI):

286053

UNIPROT ID:

Q96MF7

Full Name:

non-SMC element 2, MMS21 homolog  
(S. cerevisiae)

Calculated MW:

247 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

242862G6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01897-1, NSMCE2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85188-2-PBS. Detection antibody: 85188-1-PBS. Standard: Ag4545. Range: 1.563-200 ng/mL.