

For Research Use Only

# Caveolin-3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85190-1-PBS



## Basic Information

Catalog Number:

85190-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28009

GenBank Accession Number:

BC069368

GeneID (NCBI):

859

UNIPROT ID:

P56539

Full Name:

caveolin 3

Calculated MW:

151 aa, 17 kDa

Observed MW:

20 kDa

Purification Method:

Protein A purification

CloneNo.:

242897B2

## Applications

Tested Applications:

WB, Cytometric bead array, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85190-1-RR (Caveolin-3 antibody) at dilution of 1:30000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85190-1-PBS in a different storage buffer formulation.



Cytometric bead array standard curve of MP01900-2, Caveolin-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85190-2-PBS. Detection antibody: 85190-1-PBS. Standard: Ag28009. Range: 0.781-100 ng/mL