

For Research Use Only

p150 glued Recombinant antibody, PBS Only (Capture)

Catalog Number: 85204-1-PBS



Basic Information

Catalog Number:

85204-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004082

GeneID (NCBI):

1639

UNIPROT ID:

Q14203

Full Name:

dynactin 1 (p150, glued homolog, Drosophila)

Calculated MW:

142 kDa

Purification Method:

Protein A purification

CloneNo.:

242675A3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01959-1, Human DCTN1/p150 glued Recombinant Matched Antibody Pair - PBS only. 85204-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY01243. 85204-2-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.

The mean DCTN1/p150 glued concentration was determined to be 44,659.4 pg/mL in Jurkat cell extract based on a 3.1 mg/mL extract load.