

For Research Use Only

XDH Recombinant antibody, PBS Only (Detector)

Catalog Number: 85206-2-PBS



Basic Information

Catalog Number:

85206-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000379

GeneID (NCBI):

7498

UNIPROT ID:

P47989

Full Name:

xanthine dehydrogenase

Calculated MW:

146 kDa

Purification Method:

Protein A purification

CloneNo.:

242724F1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01961-1, Human XDH Recombinant Matched Antibody Pair - PBS only. 85206-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY01209. 85206-2-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



The mean XDH concentration was determined to be 36.55 ng/mL in mouse liver tissue extract based on a 3.7 mg/mL extract load.