

For Research Use Only

Mouse Cd72 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:85208-1-PBS



Basic Information

Catalog Number: 85208-1-PBS	GenBank Accession Number: NM_007654.2	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 12517	CloneNo.: 242846B11
Isotype: IgG	UNIPROT ID: P21855	
Immunogen Catalog Number: EG2708	Full Name: CD72 antigen	
	Calculated MW: 40 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD72 concentration was determined to be 24.3 ng/mL in RAW264.7 cell extract based on a 1.5 mg/mL extract load. Mouse brain tissue extract were served as a negative control.



Cytometric bead array standard curve of MP01912-2, MOUSE CD72 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85208-1-PBS. Detection antibody: 85208-4-PBS. Standard: Eg2708. Range: 0.781-100 ng/mL