

For Research Use Only

# Mouse Cd72 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:85208-2-PBS



## Basic Information

Catalog Number:

85208-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2708

GenBank Accession Number:

NM\_007654.2

GeneID (NCBI):

12517

UNIPROT ID:

P21855

Full Name:

CD72 antigen

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

242846D10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean CD72 concentration was determined to be 11.7 ng/mL in RAW264.7 cell extract based on a 1.5 mg/mL extract load. Mouse brain tissue extract were served as a negative control.



Cytometric bead array standard curve of MP01912-1, MOUSE CD72 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85208-3-PBS. Detection antibody: 85208-2-PBS. Standard: Eg2708. Range: 0.781-100 ng/mL