

For Research Use Only

Mouse Cd72 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:85208-3-PBS



Basic Information

Catalog Number:

85208-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2708

GenBank Accession Number:

NM_007654.2

GeneID (NCBI):

12517

UNIPROT ID:

P21855

Full Name:

CD72 antigen

Calculated MW:

40 kDa

Purification Method:

Protein A purification

CloneNo.:

242846E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD72 concentration was determined to be 11.7 ng/mL in RAW264.7 cell extract based on a 1.5 mg/mL extract load. Mouse brain tissue extract were served as a negative control.



Cytometric bead array standard curve of MP01912-1, MOUSE CD72 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85208-3-PBS. Detection antibody: 85208-2-PBS. Standard: Eg2708. Range: 0.781-100 ng/mL