

For Research Use Only

CD40L/CD154 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85209-1-PBS



Basic Information

Catalog Number:

85209-1-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3501

GenBank Accession Number:

BC071754

GeneID (NCBI):

959

UNIPROT ID:

P29965

Full Name:

CD40 ligand

Calculated MW:

261 aa, 29 kDa

Purification Method:

Protein A purification

CloneNo.:

242847B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01907-2, CD40L/CD154 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85209-2-PBS. Detection antibody: 85209-1-PBS. Standard: Eg3501. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 85209-1-PBS against Human CD40L/CD154 was performed. The affinity constant is 1.07 nM.