

For Research Use Only

# KCTD5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85210-1-PBS



## Basic Information

Catalog Number:

85210-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7915

GenBank Accession Number:

BC007314

GeneID (NCBI):

54442

UNIPROT ID:

Q9NXV2

Full Name:

potassium channel tetramerisation  
domain containing 5

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

242864A10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01913-2, KCTD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85210-1-PBS. Detection antibody: 85210-2-PBS. Standard: Ag7915. Range: 0.781-100 ng/mL.