

For Research Use Only

KCTD5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85210-2-PBS



Basic Information

Catalog Number:

85210-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7915

GenBank Accession Number:

BC007314

GeneID (NCBI):

54442

UNIPROT ID:

Q9NXV2

Full Name:

potassium channel tetramerisation
domain containing 5

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

242864A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01913-2, KCTD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85210-1-PBS. Detection antibody: 85210-2-PBS. Standard: Ag7915. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01913-1, KCTD5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85210-3-PBS. Detection antibody: 85210-2-PBS. Standard: Ag7915. Range: 0.781-100 ng/mL