

For Research Use Only

Rat MMP8 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85211-1-PBS



Basic Information

Catalog Number:

85211-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_022221.1

GeneID (NCBI):

63849

UNIPROT ID:

O88766

Full Name:

matrix metalloproteinase 8

Calculated MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

242473B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

Background Information

MMP8(Matrix metalloproteinase-8) an enzyme that degrades fibrillar collagens imparting strength to the fetal membranes, is expressed by leukocytes and chorionic cytotrophoblast cells(PMID:15367487).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01911-2, RAT MMP-8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85211-1-PBS. Detection antibody: 85211-3-PBS. Standard: Eg3278. Range: 0.469-60 ng/mL.



Cytometric bead array standard curve of MP01911-1, RAT MMP-8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85211-1-PBS. Detection antibody: 85211-2-PBS. Standard: Eg3278. Range: 0.469-60 ng/mL.