

For Research Use Only

ATRIP Recombinant antibody, PBS Only (Detector)

Catalog Number: 85213-2-PBS



Basic Information

Catalog Number:

85213-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1869

GenBank Accession Number:

BC014153

GeneID (NCBI):

84126

UNIPROT ID:

Q8WXE1

Full Name:

ATR interacting protein

Calculated MW:

791 aa, 86 kDa

Purification Method:

Protein A purification

CloneNo.:

242756C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01905-2, ATRIP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85213-3-PBS. Detection antibody: 85213-2-PBS. Standard: Ag1869. Range: 0.625-80 ng/mL