

For Research Use Only

# SLC39A8 Recombinant antibody, PBS Only (Detector)

Catalog Number:85215-4-PBS



## Basic Information

Catalog Number:

85215-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14815

GenBank Accession Number:

BC012125

GeneID (NCBI):

64116

UNIPROT ID:

Q9C0K1

Full Name:

solute carrier family 39 (zinc transporter), member 8

Calculated MW:

460 aa, 50 kDa

Purification Method:

Protein A purification

CloneNo.:

242773H8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01908-2, ZIP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85215-3-PBS. Detection antibody: 85215-4-PBS. Standard: Ag14815. Range: 1.563-100 ng/mL.