

For Research Use Only

DLL4 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85216-3-PBS



Basic Information

Catalog Number:

85216-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16140

GenBank Accession Number:

BC106950

GeneID (NCBI):

54567

UNIPROT ID:

Q9NR61

Full Name:

delta-like 4 (Drosophila)

Calculated MW:

685 aa, 75 kDa

Purification Method:

Protein A purification

CloneNo.:

242833G3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01916-1, DLL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85216-1-PBS. Detection antibody: 85216-3-PBS. Standard: Ag16140. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01916-2, DLL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85216-2-PBS. Detection antibody: 85216-3-PBS. Standard: Ag16140. Range: 0.781-100 ng/mL.