

For Research Use Only

# DFFA,DFF45 Recombinant antibody, PBS Only (Detector)

Catalog Number:85220-2-PBS



## Basic Information

Catalog Number:

85220-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0244

GenBank Accession Number:

BC000037

GeneID (NCBI):

1676

UNIPROT ID:

O00273

Full Name:

DNA fragmentation factor, 45kDa,  
alpha polypeptide

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

242915E12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01914-2, DFFA/DFF45 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85220-3-PBS. Detection antibody: 85220-2-PBS. Standard: Ag0244. Range: 0.781-100 ng/mL.