

For Research Use Only

DFFA,DFF45 Recombinant antibody, PBS Only (Capture)

Catalog Number:85220-3-PBS



Basic Information

Catalog Number:

85220-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0244

GenBank Accession Number:

BC000037

GeneID (NCBI):

1676

UNIPROT ID:

O00273

Full Name:

DNA fragmentation factor, 45kDa,
alpha polypeptide

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

242915G9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01914-2, DFFA/DFF45 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85220-3-PBS. Detection antibody: 85220-2-PBS. Standard: Ag0244. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01914-1, DFFA/DFF45 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85220-3-PBS. Detection antibody: 85220-1-PBS. Standard: Ag0244. Range: 0.781-100 ng/mL.