

For Research Use Only

SCN9A Recombinant antibody, PBS Only (Capture)

Catalog Number: 85279-2-PBS



Basic Information

Catalog Number:

85279-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001365536.1

GeneID (NCBI):

6335

UNIPROT ID:

Q15858-1

Full Name:

sodium channel, voltage-gated, type IX, alpha subunit

Calculated MW:

226kDa

Purification Method:

Protein A purification

CloneNo.:

242997F10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02070-2, SCN9A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85279-2-PBS. Detection antibody: 85279-3-PBS. Standard: SY02473. Range: 1.563-200 ng/mL.