

For Research Use Only

FGFBasic-TS Recombinant antibody, PBS Only (Capture)

Catalog Number: 85284-3-PBS



Basic Information

Catalog Number:

85284-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

HZ-1285

GenBank Accession Number:

GeneID (NCBI):

2247

UNIPROT ID:

P09038

Full Name:

fibroblast growth factor 2 (basic)

Purification Method:

Protein A purification

CloneNo.:

250430E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FGF2 concentration was determined to be 113.2 ng/mL in HeLa cell extract based on a 1.3 mg/mL extract load.



Cytometric bead array standard curve of MP02231-1, FGF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85284-3-PBS. Detection antibody: 85284-2-PBS. Standard: HZ-1285. Range: 0.781-50 ng/mL.