

For Research Use Only

NID2 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:85299-3-PBS



Basic Information

Catalog Number: 85299-3-PBS	GenBank Accession Number: BC035608	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 22795	CloneNo.: 25235088
Isotype: IgG	UNIPROT ID: Q14112	
Immunogen Catalog Number: AG4643	Full Name: nidogen 2 (osteonidogen)	
	Calculated MW: 1375 aa, 151 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Nidogen 2 concentration was determined to be 49.3 ng/mL in human placenta cell extract based on a 3.2 mg/mL extract load, HT-29 cell extract were served as a negative control.



Cytometric bead array standard curve of MP02894-1, Nidogen 2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85299-3-PBS. Detection antibody: 85299-2-PBS. Standard: Ag4643. Range: 1.25-80 ng/mL.