

For Research Use Only

MYH7 Recombinant antibody

Catalog Number: 85322-1-RR



Basic Information

Catalog Number:

85322-1-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000257

GeneID (NCBI):

4625

UNIPROT ID:

P12883

Full Name:

myosin, heavy chain 7, cardiac muscle, beta

Calculated MW:

223 kDa

Observed MW:

200-220 kDa

Purification Method:

Protein A purification

CloneNo.:

242250F6

Recommended Dilutions:

WB 1:2000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : mouse heart tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse heart tissue were subjected to SDS PAGE followed by western blot with 85322-1-RR (MYH7 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.