

For Research Use Only

DTX3L Recombinant antibody, PBS Only

Catalog Number: 85340-2-PBS



Basic Information

Catalog Number:

85340-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2568

GenBank Accession Number:

BC060509

GeneID (NCBI):

151636

UNIPROT ID:

Q8TDB6

Full Name:

deltex 3-like (Drosophila)

Calculated MW:

740 aa, 84 kDa

Observed MW:

84 kDa

Purification Method:

Protein A purification

CloneNo.:

242627C8

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85340-2-RR (DTX3L antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85340-2-PBS in a different storage buffer formulation.