

For Research Use Only

TRIM33 Recombinant antibody, PBS Only

Catalog Number: 85364-3-PBS



Basic Information

Catalog Number:

85364-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_015906

GeneID (NCBI):

51592

UNIPROT ID:

Q9UPN9

Full Name:

tripartite motif-containing 33

Calculated MW:

122 kDa

Observed MW:

140-150 kDa

Purification Method:

Protein A purification

CloneNo.:

242728H2

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

human

Background Information

TRIM33, also named KIAA1113, RFG7, and TIF1G, belongs to the TRIM/RBCC family. It acts as an E3 ubiquitin-protein ligase. TRIM33 promotes SMAD4 ubiquitination, nuclear exclusion, and degradation via the ubiquitin-proteasome pathway. It may act as a transcriptional repressor. It inhibits the transcriptional response to TGF-beta/BMP signaling cascade. TRIM33 plays a role in the control of cell proliferation. Its association with SMAD2 and SMAD3, TRIM33 stimulates erythroid differentiation of hematopoietic stem/progenitor. Defects in TRIM33 are a cause of thyroid papillary carcinoma (TPC). This antibody is specific to TRIM33. The calculated molecular weight of TRIM33 is 122kDa, but the molecular weight detected is about 140-150 kDa the same as the result from the paper(PMID:25639486;26592194).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85364-3-RR (TRIM33 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85364-3-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using TRIM33 antibody (85364-3-RR, Clone: 242728H2) at dilution of 1:500 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 85364-3-PBS in a different storage buffer formulation.



COLO 320 cells were subjected to SDS PAGE followed by western blot with 85364-3-RR (TRIM33 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85364-3-PBS in a different storage buffer formulation.