

For Research Use Only

MYH3 Recombinant monoclonal antibody

Catalog Number: 85369-1-RR



Basic Information

Catalog Number:

85369-1-RR

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002470

GeneID (NCBI):

4621

UNIPROT ID:

P11055

Full Name:

myosin, heavy chain 3, skeletal muscle, embryonic

Calculated MW:

234 kDa

Observed MW:

220-230 kDa

Purification Method:

Protein A purification

CloneNo.:

242722C1

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : mouse embryo tissue, mouse eye tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85369-1-RR (MYH3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.