

For Research Use Only

MAP4K4 Recombinant monoclonal antibody

Catalog Number: 85387-1-RR



Basic Information

Catalog Number:

85387-1-RR

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_004834

GeneID (NCBI):

9448

UNIPROT ID:

O95819

Full Name:

mitogen-activated protein kinase
kinase kinase kinase 4

Calculated MW:

151 kDa

Purification Method:

Protein A purification

CloneNo.:

242721G1

Recommended Dilutions:

IF/ICC: 1:200-1:800

FC (Intra): 0.25 ug per 10⁶ cells in a
100 µl suspension

Applications

Tested Applications:

IF/ICC, FC (Intra), ELISA

Species Specificity:

human

Positive Controls:

IF/ICC : HeLa cells,

FC (Intra) : HepG2 cells,

Background Information

MAP4K4, also named as RAB8IPL1, GLK and MEKKK4, belongs to the protein kinase superfamily, STE Ser/Thr protein kinase family and STE20 subfamily. It may play a role in the response to environmental stress. MAP4K4 appears to act upstream of the JUN N-terminal pathway.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using MAP4K4 antibody (85387-1-RR, Clone: 242721G1) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).



1x10⁶ HepG2 cells were intracellularly stained with 0.25 ug MAP4K4 Recombinant monoclonal antibody (85387-1-RR, Clone:242721G1) and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2)(red), or 0.25 ug Rabbit IgG Isotype Control RecAb (98136-1-RR, Clone: 240953C9) (blue). Cells were fixed and permeabilized with Intracellular Flow Cytometry Fixation & Permeabilization Buffer Kit (PF00019).