

For Research Use Only

# RORB Recombinant antibody, PBS Only (Capture)

Catalog Number: 85404-3-PBS



## Basic Information

Catalog Number:

85404-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11910

GenBank Accession Number:

BC051830

GeneID (NCBI):

6096

UNIPROT ID:

Q92753

Full Name:

RAR-related orphan receptor B

Calculated MW:

459 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

242832G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01924-2, RORB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85404-3-PBS. Detection antibody: 85404-1-PBS. Standard: Ag11910. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP01924-1, RORB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85404-3-PBS. Detection antibody: 85404-2-PBS. Standard: Ag11910. Range: 1.563-200 ng/mL.