

For Research Use Only

SLC5A7 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85405-2-PBS



Basic Information

Catalog Number:

85405-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16136

GenBank Accession Number:

BC111525

GeneID (NCBI):

60482

UNIPROT ID:

Q9GZV3

Full Name:

solute carrier family 5 (choline transporter), member 7

Calculated MW:

580 aa, 63 kDa

Purification Method:

Protein A purification

CloneNo.:

242904G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01923-2, CHT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85405-2-PBS. Detection antibody: 85405-3-PBS. Standard: Ag16136. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01923-1, CHT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85405-2-PBS. Detection antibody: 85405-1-PBS. Standard: Ag16136. Range: 0.781-100 ng/mL.