

For Research Use Only

# PTGES2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85419-1-PBS



## Basic Information

Catalog Number:

85419-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1329

GenBank Accession Number:

BC011613

GeneID (NCBI):

80142

UNIPROT ID:

Q9H7Z7

Full Name:

prostaglandin E synthase 2

Calculated MW:

42 kDa

Purification Method:

Protein A purification

CloneNo.:

242930A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01948-2, PTGES2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85419-2-PBS. Detection antibody: 85419-1-PBS. Standard: Ag1329. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01948-1, PTGES2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85419-1-PBS. Detection antibody: 85419-3-PBS. Standard: Ag1329. Range: 0.781-100 ng/mL.