

For Research Use Only

FGFR4 Recombinant monoclonal antibody, PBS Only

Catalog Number: 85420-3-PBS



Basic Information

Catalog Number: 85420-3-PBS	GenBank Accession Number: BC011847	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 2264	CloneNo.: 242713H8
Isotype: IgG	UNIPROT ID: P22455	
Immunogen Catalog Number: EG3510	Full Name: fibroblast growth factor receptor 4	
	Calculated MW: 88 kDa	
	Observed MW: 95 kDa, 110 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

Fibroblast growth factor receptor 4 (FGFR4) is a member of a highly conserved tyrosine kinase family, along with FGFR1-3. This family consists of an intracellular tyrosine kinase domain, a single transmembrane domain, and extracellular ligand-binding domains (PMID: 32492514). FGFR4 is the predominant FGFR isoform present in human hepatocytes. FGFR4 has been proposed to play a role in the observed induction of hepatocyte proliferation and carcinogenesis by FGF19; however, contradicting evidence proposing a protective role for FGFR4 in suppressing hepatoma progression has also been proposed (PMID: 20018895). While the role of FGFR4 in cancer remains to be fully elucidated, several findings suggest that this receptor may be an important player in Hepatocellular carcinoma (HCC) development and/or progression (PMID: 10336501).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85420-3-RR (FGFR4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85420-3-PBS in a different storage buffer formulation.