

For Research Use Only

KRT6B Recombinant antibody, PBS Only (Capture)

Catalog Number: 85423-3-PBS



Basic Information

Catalog Number:

85423-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11288

GenBank Accession Number:

BC034535

GeneID (NCBI):

3854

UNIPROT ID:

P04259

Full Name:

keratin 6B

Calculated MW:

564 aa, 60 kDa

Purification Method:

Protein A purification

CloneNo.:

242874G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01943-1, Cytokeratin 6B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85423-3-PBS. Detection antibody: 85423-1-PBS. Standard: Ag11288. Range: 0.781-100 ng/mL.