

For Research Use Only

FUT11 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85427-3-PBS



Basic Information

Catalog Number:

85427-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10815

GenBank Accession Number:

BC036037

GeneID (NCBI):

170384

UNIPROT ID:

Q495W5

Full Name:

fucosyltransferase 11 (alpha (1,3)

fucosyltransferase)

Calculated MW:

492 aa, 56 kDa

Purification Method:

Protein A purification

CloneNo.:

242869H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01944-1, FUT11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85427-1-PBS. Detection antibody: 85427-3-PBS. Standard: Ag10815. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01944-2, FUT11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85427-2-PBS. Detection antibody: 85427-3-PBS. Standard: Ag10815. Range: 0.625-80 ng/mL