

For Research Use Only

# ILK Recombinant antibody, PBS Only (Capture)

Catalog Number: 85428-3-PBS



## Basic Information

Catalog Number: 85428-3-PBS	GenBank Accession Number: BC001554	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 3611	CloneNo.: 242982F5
Source: Rabbit	UNIPROT ID: Q13418	
Isotype: IgG	Full Name: integrin-linked kinase	
Immunogen Catalog Number: AG3570	Calculated MW: 451 aa, 59 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01945-2, ILK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85428-3-PBS. Detection antibody: 85428-4-PBS. Standard: Ag3570. Range: 0.781-100 ng/mL.