

For Research Use Only

PEX1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85429-3-PBS



Basic Information

Catalog Number:

85429-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4623

GenBank Accession Number:

BC035575

GeneID (NCBI):

5189

UNIPROT ID:

O43933

Full Name:

peroxisomal biogenesis factor 1

Calculated MW:

1283 aa, 143 kDa

Purification Method:

Protein A purification

CloneNo.:

242994G11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01946-1, PEX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85429-1-PBS. Detection antibody: 85429-3-PBS. Standard: Ag4623. Range: 0.781-100 ng/mL.