

For Research Use Only

PEX1 Recombinant antibody, PBS Only

Catalog Number: 85429-4-PBS



Basic Information

Catalog Number:

85429-4-PBS

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4623

GenBank Accession Number:

BC035575

GeneID (NCBI):

5189

UNIPROT ID:

O43933

Full Name:

peroxisomal biogenesis factor 1

Calculated MW:

1283 aa, 143 kDa

Observed MW:

143 kDa

Purification Method:

Protein A purification

CloneNo.:

242994F11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

Background Information

Peroxin (PEX) proteins generate and maintain peroxisomes. The peroxin-1 (PEX1) ATPase facilitates the recycling of the peroxisome matrix protein receptor PEX5 and is the most commonly affected peroxin in human peroxisome biogenesis disorders (PMID: 28600347). PEX1 and PEX6 are the only members of the AAA family (for ATPases associated with diverse cellular) activities implicated in peroxisome biogenesis and are closely related to p97, which functions in endoplasmic reticulum-associated protein degradation to retrotranslocate endoplasmic reticulum proteins to the cytosol (PMID: 9695811; 11740563).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85429-4-RR (PEX1 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85429-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85429-4-RR against Human PEX1 were performed. The affinity constant is below 1 pM.