

For Research Use Only

# CABC1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85431-3-PBS



## Basic Information

Catalog Number:

85431-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7612

GenBank Accession Number:

BC005171

GeneID (NCBI):

56997

UNIPROT ID:

Q8NI60

Full Name:

chaperone, ABC1 activity of bc1 complex homolog (S. pombe)

Calculated MW:

72 kDa

Purification Method:

Protein A purification

CloneNo.:

242942F10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01952-1, CABCl Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85431-3-PBS. Detection antibody: 85431-1-PBS. Standard: Ag7612. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01952-2, CABCl Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85431-3-PBS. Detection antibody: 85431-2-PBS. Standard: Ag7612. Range: 0.781-100 ng/mL