

For Research Use Only

CTR9 Recombinant antibody, PBS Only

Catalog Number: 85432-4-PBS



Basic Information

Catalog Number:

85432-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15690

GenBank Accession Number:

BC058914

GeneID (NCBI):

9646

UNIPROT ID:

Q6PD62

Full Name:

Ctr9, Paf1/RNA polymerase II complex component, homolog (S. cerevisiae)

Calculated MW:

1173 aa, 134 kDa

Observed MW:

140 kDa

Purification Method:

Protein A purification

CloneNo.:

242949B10

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

CTR9, also named as KIAA0155 and SH2BP1, contains 16 TPR repeats. It is widely expressed.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85432-4-RR (CTR9 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85432-4-PBS in a different storage buffer formulation.