

For Research Use Only

# CTR9 Recombinant antibody

Catalog Number: 85432-4-RR



## Basic Information

Catalog Number:

85432-4-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15690

GenBank Accession Number:

BC058914

GeneID (NCBI):

9646

UNIPROT ID:

Q6PD62

Full Name:

Ctr9, Paf1/RNA polymerase II complex component, homolog (S. cerevisiae)

Calculated MW:

1173 aa, 134 kDa

Observed MW:

140 kDa

Purification Method:

Protein A purification

CloneNo.:

242949B10

Recommended Dilutions:

WB 1:5000-1:50000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HEK-293 cells, L02 cells, mouse lung tissue, rat testis tissue

## Background Information

CTR9, also named as KIAA0155 and SH2BP1, contains 16 TPR repeats. It is widely expressed.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85432-4-RR (CTR9 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.