

For Research Use Only

ABCA1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85434-1-PBS



Basic Information

Catalog Number:

85434-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24055

GenBank Accession Number:

BC141816

GeneID (NCBI):

19

UNIPROT ID:

O95477

Full Name:

ATP-binding cassette, sub-family A
(ABC1), member 1

Calculated MW:

2261 aa, 254 kDa

Purification Method:

Protein A purification

CloneNo.:

242872C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01937-1, ABCA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85434-2-PBS. Detection antibody: 85434-1-PBS. Standard: Ag24055. Range: 1.563-200 ng/mL.