

For Research Use Only

# Mouse Cdkn1a Recombinant antibody, PBS Only (Capture)

Catalog Number: 85435-1-PBS



## Basic Information

Catalog Number:

85435-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28249

GenBank Accession Number:

NM\_001111099

GeneID (NCBI):

12575

UNIPROT ID:

P39689

Full Name:

cyclin-dependent kinase inhibitor 1A  
(P21)

Calculated MW:

18 kDa

Purification Method:

Protein A purification

CloneNo.:

242478E9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01951-1, MOUSE p21 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85435-1-PBS. Detection antibody: 85435-2-PBS. Standard:Ag28249. Range: 0.781-100 ng/mL.