

For Research Use Only

CD158G/KIR2DS5 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85436-1-PBS



Basic Information

Catalog Number:

85436-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_014513.3

GeneID (NCBI):

3810

UNIPROT ID:

Q14953

Full Name:

killer cell immunoglobulin-like receptor, two domains, short cytoplasmic tail, 5

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

242877A11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01954-1, CD158G/KIR2DS5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85436-1-PBS. Detection antibody: 85436-2-PBS. Standard: Eg3719. Range: 0.781-100 ng/mL.