

For Research Use Only

GHR Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85437-1-PBS



Basic Information

Catalog Number:

85437-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3027

GenBank Accession Number:

NM_000163.5

GeneID (NCBI):

2690

UNIPROT ID:

P10912-1

Full Name:

GH binding protein

Calculated MW:

72 kDa

Purification Method:

Protein A purification

CloneNo.:

242896B8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01942-2, GHR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85437-3-PBS. Detection antibody: 85437-1-PBS. Standard: Eg3027. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01942-1, GHR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85437-1-PBS. Detection antibody: 85437-2-PBS. Standard: Eg3027. Range: 0.313-40 ng/mL