

For Research Use Only

GHR Recombinant antibody, PBS Only (Detector)

Catalog Number: 85437-4-PBS



Basic Information

Catalog Number:

85437-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3027

GenBank Accession Number:

NM_000163.5

GeneID (NCBI):

2690

UNIPROT ID:

P10912-1

Full Name:

GH binding protein

Calculated MW:

72 kDa

Purification Method:

Protein A purification

CloneNo.:

242896F3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01942-3, Human GHR Recombinant Matched Antibody Pair - PBS only. 85437-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3027. 85437-4-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



Serum of six individual healthy human donors was measured. The GHR concentration of detected samples was determined to be 18.0 ng/mL with a range of 6.6-25.6 ng/mL.