

For Research Use Only

CHGB Recombinant antibody, PBS Only (Capture)

Catalog Number: 85438-1-PBS



Basic Information

Catalog Number: 85438-1-PBS	GenBank Accession Number: BC000375	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 1114	CloneNo.: 242952B6
Source: Rabbit	UNIPROT ID: P05060	
Isotype: IgG	Full Name: chromogranin B (secretogranin 1)	
Immunogen Catalog Number: AG6833	Calculated MW: 78 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01936-2, Chromogranin B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85438-1-PBS. Detection antibody: 85438-3-PBS. Standard: Ag6833. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01936-1, Chromogranin B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85438-1-PBS. Detection antibody: 85438-2-PBS. Standard: Ag6833. Range: 0.781-100 ng/mL.