

For Research Use Only

CD96 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85439-5-PBS



Basic Information

Catalog Number:

85439-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3520

GenBank Accession Number:

NM_198196

GeneID (NCBI):

10225

UNIPROT ID:

P40200

Full Name:

CD96 molecule

Calculated MW:

402 aa, 46 kDa

Purification Method:

Protein A purification

CloneNo.:

242964F9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01941-2, CD96 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85439-3-PBS. Detection antibody: 85439-5-PBS. Standard: Eg3520. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 85439-5-PBS against Human CD96 was performed. The affinity constant is 4.25 nM.