

For Research Use Only

# MSR1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85441-2-PBS



## Basic Information

Catalog Number:

85441-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC063878

GeneID (NCBI):

4481

UNIPROT ID:

P21757

Full Name:

macrophage scavenger receptor 1

Calculated MW:

451 aa, 50 kDa

Purification Method:

Protein A purification

CloneNo.:

242710E3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01935-2, CD204 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85441-2-PBS. Detection antibody: 85441-1-PBS. Standard: Eg2566. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP01935-1, CD204 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85441-2-PBS. Detection antibody: 85441-3-PBS. Standard: Eg2566. Range: 1.563-200 ng/mL.