

For Research Use Only

# HPGDS Recombinant antibody, PBS Only (Capture)

Catalog Number: 85446-3-PBS



## Basic Information

Catalog Number:

85446-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18290

GenBank Accession Number:

BC020734

GeneID (NCBI):

27306

UNIPROT ID:

O60760

Full Name:

prostaglandin D2 synthase,  
hematopoietic

Calculated MW:

199 aa, 23 kDa

Purification Method:

Protein A purification

CloneNo.:

242980D6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP02001-2, HPGDS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85446-3-PBS. Detection antibody: 85446-2-PBS. Standard: Ag18290. Range: 0.781-50 ng/mL