

For Research Use Only

SLC32A1/VGAT Recombinant antibody, PBS Only (Capture)

Catalog Number: 85447-2-PBS



Basic Information

Catalog Number:

85447-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5843

GenBank Accession Number:

BC053582

GeneID (NCBI):

140679

UNIPROT ID:

Q9H598

Full Name:

solute carrier family 32 (GABA vesicular transporter), member 1

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

242990F8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01938-2, VGAT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85447-2-PBS. Detection antibody: 85447-3-PBS. Standard: Ag5843. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01938-1, VGAT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85447-2-PBS. Detection antibody: 85447-1-PBS. Standard: Ag5843. Range: 0.781-100 ng/mL.